---

# EIAgen

## HBsAg Kit

---

**REF** 071011


$\Sigma$  96

**REF** 071012

$\Sigma$  192

**REF** 071015

$\Sigma$  480

**IVD**

This package insert must be read carefully before product use.

Package insert instructions must be carefully followed.

Reliability of assay results cannot be guaranteed if there are any deviations from the instructions in this package.

**Manufacturer:**

**Adaltis S.r.l**

Via Durini, 27

20122 Milano (Italy)

Tel. +39-0774-5791 - Fax +39-0774-353085

[www.adaltis.net](http://www.adaltis.net)

**en**

**SYMBOLS USED ON LABELS**

<b>English</b> <b>EN</b>							
	In Vitro Diagnostic Medical Device	Catalogue Number	Lot Number	Attention, See Instructions For Use	Temperature Limitation	Use By	Number of Test
							
	Manufacturer	Keep away from Sunlight	Date of Manufacture	Biological Risk	Microplate	Positive Control	Negative Control
							
	Conjugate	Assay Diluent	Substrate A (Urea Peroxide)	Substrate B (TMB)	Stop Solution (0,5 M H <sub>2</sub> SO <sub>4</sub> )	Wash Buffer Conc. (20x)	

**A. INTENDED USE**

EIAgen HBsAg Kit is an enzyme-linked immunosorbent assay (ELISA) for qualitative detection of HBsAg in human serum or plasma (EDTA, sodium citrate or heparin). It is intended for screening of blood donors and for diagnosing of patients related to infection with hepatitis B virus.

For "in vitro" diagnostic use only.

**B. INTRODUCTION**

Hepatitis B virus (HBV) is an enveloped, double-stranded DNA virus belonging to the Hepadnaviridae family and is recognized as the major cause of blood transmitted hepatitis together with hepatitis C virus (HCV). Infection with HBV induces a spectrum of clinical manifestations ranging from mild, inapparent disease to fulminant hepatitis, severe chronic liver diseases, which in some cases can lead to cirrhosis and carcinoma of the liver. Classification of a hepatitis B infection requires the identification of several serological markers expressed during three phases (incubation, acute and convalescent) of the infection. Now several diagnostic test are used for screening, clinical diagnosis and management of the disease. Hepatitis B surface antigen or HBsAg, previously described as Australia antigen, is the most important protein of the envelope of Hepatitis B Virus. The surface antigen contains the determinant "a", common to all known viral subtypes and immunologically distinguished in two distinct subgroups (ay and ad). HBV has 10 major serotypes and four HBsAg subtypes have been recognized (adw, ady, ayw, and ayr). HBsAg can be detected 2 to 4 weeks before the ALT levels become abnormal and 3 to 5 weeks before symptoms develop. The serological detection of HBsAg is a powerful method for the diagnosis and prevention of HBV infection and ELISA has become an extensively used analytical system for screening of blood donors and clinical diagnosis of HBV in infected individuals.

**C. PRINCIPLE OF THE TEST**

For detection of HBsAg, EIAgen HBsAg Kit uses antibody "sandwich" ELISA method in which, polystyrene microwell strips are pre-coated with monoclonal antibodies specific to HBsAg. Patient's serum or plasma sample is added to the microwells. During incubation, the specific immunocomplex formed in case of presence of HBsAg in the sample, is captured on the solid phase. Then the second antibody conjugated the enzyme horseradish peroxidase (the HRP-Conjugate) directed against a different epitope of HBsAg is added into the wells. During the second incubation step, these HRP-conjugated antibodies will be bound to any anti-HBs-HBsAg complexes previously formed during the first incubation, and the unbound HRP-conjugate is then removed by washing. Chromogen solutions containing tetramethyl-benzidine (TMB) and urea peroxide are added to the wells. In presence of the antibody-antigen-antibody(HRP) "sandwich" immunocomplex, the colorless Chromogens are hydrolyzed by the bound HRP-conjugate to a blue-colored product. The blue color turns yellow after stopping the reaction with sulfuric acid. The amount of color intensity can be measured and it is proportional to the amount of antigen captured in the

wells, and to its amount in the sample respectively. Wells containing samples negative for HBsAg remain colorless.

**D. COMPONENTS**

The kit contains reagents for 96 tests (code 071011), 192 tests (code 071012), or 480 tests (code 071015).

<b>Microplate</b>	1
<b>Negative Control</b>	1 x 1.5 mL/vial
<b>Positive Control</b>	1 x 1.5 mL/vial
<b>Conjugate</b>	1 x 6 mL/vial
<b>Assay Diluent</b>	1 x 5 mL/vial
<b>Substrate Solution A (Urea Peroxide)</b>	1 x 6 mL/vial
<b>Substrate Solution B (TMB)</b>	1 x 6 mL/vial
<b>Stop Solution (H<sub>2</sub>SO<sub>4</sub> - 0.5M)</b>	1 x 6 mL/vial
<b>Wash Buffer Concentrate 20X</b>	1 x 30 mL/vial
<b>Plate Sealing Foils</b>	2
<b>Plastic Sealable Bag</b>	1
<b>Number of tests</b>	96
<b>Code</b>	071011

<b>Microplate</b>	2
<b>Negative Control</b>	1 x 1.5 mL/vial
<b>Positive Control</b>	1 x 1.5 mL/vial
<b>Conjugate</b>	2 x 6 mL/vial
<b>Assay Diluent</b>	2 x 5 mL/vial
<b>Substrate Solution A (Urea Peroxide)</b>	2 x 6 mL/vial
<b>Substrate Solution B (TMB)</b>	2 x 6 mL/vial
<b>Stop Solution (H<sub>2</sub>SO<sub>4</sub> - 0.5M)</b>	2 x 6 mL/vial
<b>Wash Buffer Concentrate 20X</b>	2 x 30 mL/vial
<b>Plate Sealing Foils</b>	3
<b>Plastic Sealable Bag</b>	2
<b>Number of tests</b>	192
<b>Code</b>	071012

<b>Microplate</b>	5
<b>Negative Control</b>	4 x 1.5 mL/vial
<b>Positive Control</b>	4 x 1.5 mL/vial
<b>Conjugate</b>	5 x 6 mL/vial
<b>Assay Diluent</b>	4 x 5 mL/vial
<b>Substrate Solution A (Urea Peroxide)</b>	5 x 6 mL/vial
<b>Substrate Solution B (TMB)</b>	5 x 6 mL/vial
<b>Stop Solution (H<sub>2</sub>SO<sub>4</sub> - 0.5M)</b>	1 x 30 mL/vial
<b>Wash Buffer Concentrate 20X</b>	3 x 50 mL/vial
<b>Plate Sealing Foils</b>	6
<b>Plastic Sealable Bag</b>	5
<b>Number of tests</b>	480
<b>Code</b>	071015

**1. Microplate**

Blank microwell strips fixed on white strip holder.

12 strips of 8 microwells coated with monoclonal antibodies reactive to HBsAg (anti-HBs).

Plates are sealed into a aluminium pouch with desiccant. The microwell strips can be broken to be used separately. Place unused wells or strips in the provided plastic sealable storage bag together with the desiccant and return to 2...8°C. Once open, stable for one month at 2...8°C.

**2. Negative Control**

Yellowish liquid filled in a vial with natural screw cap.

Protein-stabilized buffer tested non-reactive for HBsAg.

Ready to use as supplied. It contains preservative.

Once open, stable for one month at 2...8°C.

### 3. Positive Control

Red-colored liquid filled in a vial with red screw cap.  
HBsAg diluted in protein-stabilized buffer.  
Ready to use as supplied. It contains preservative.  
Once open, stable for one month at 2...8°C.

**Important Note:** *The absence of viable pathogens in the Positive Control can not be fully ensured, and therefore, the reagent should be handled as potentially biohazardous, in accordance with good laboratory practices.*

### 4. Conjugate

Red-colored liquid in a white vial with green screw cap.  
Horseradish peroxidase-conjugated anti-HBs antibodies.  
Ready to use as supplied. It contains preservative.  
Once open, stable for one month at 2...8°C.

### 5. Assay Diluent

Green-colored in a vial with pink screw cap.  
Serum base, casein, and sucrose solution.  
Ready to use as supplied. It contains preservative.  
Once open, stable for one month at 2...8°C.

### 6. Substrate Solution "A"

Colorless liquid filled in a white vial with grey screw cap.  
Urea peroxide solution.  
Ready to use as supplied.  
Once open, stable for one month at 2...8°C.

### 7. Substrate Solution "B"

Colorless liquid filled in a black vial with black screw cap.  
TMB (Tetramethyl benzidine solution).  
Ready to use as supplied.  
Once open, stable for one month at 2...8°C.

**Note:** *To be stored protected from light as sensitive to strong illumination.*

### 8. Stop Solution

Colorless liquid in a white vial with red screw cap.  
Diluted sulfuric acid solution (0.5M H<sub>2</sub>SO<sub>4</sub>).  
Ready to use as supplied.  
Once open, stable for one month at 2...8°C.

**The MSDS is available upon request of laboratory personnel.**

### 9. Wash Buffer Concentrate 20x

Colorless liquid filled in a clear bottle with natural screw cap. pH 7.4, 20 × PBS  
The concentrate must be diluted **1 to 20** with distilled/deionized water before use.  
Once diluted, stable for one week at room temperature, or for two weeks when stored at 2...8°C.

## E. MATERIALS REQUIRED BUT NOT PROVIDED

1. Calibrated Micropipettes (20, 50, 100 µL) and disposable plastic tips.
2. EIA grade water (bidistilled or deionised, charcoal treated to remove oxidizing chemicals used as disinfectants).
3. Timer with 60 minute range or higher.
4. Absorbent paper tissues.
5. Calibrated ELISA microplate thermostatic incubator capable to provide a temperature of 37 ± 0.5°C.

6. Calibrated ELISA microwell reader with 450nm (reading) and possibly with 620-630nm (blanking) filters.
7. Calibrated ELISA microplate washer.
8. Vortex or similar mixing tools.

## F. WARNINGS AND PRECAUTIONS

1. The kit has to be used by skilled and properly trained technical personnel only, under the supervision of a medical doctor responsible of the laboratory.  
This package insert must be read carefully before product use.
2. When the kit is used for the screening of blood units and blood components, it has to be used in a laboratory certified and qualified by the national authority in that field (Ministry of Health or similar entity) to carry out this type of analysis.
3. All the personnel involved in performing the assay have to wear protective laboratory clothes, talc-free gloves and glasses. The use of any sharp (needles) or cutting (blades) devices should be avoided. All the personnel involved should be trained in biosafety procedures, as recommended by the Center for Disease Control, Atlanta, U.S. and reported in the National Institute of Health's publication: "Biosafety in Microbiological and Biomedical Laboratories", ed. 1984.
4. All the personnel involved in sample handling should be vaccinated for HBV and HAV, for which vaccines are available, safe and effective.
5. The laboratory environment should be controlled so as to avoid contaminants such as dust or air-borne microbial agents, when opening kit vials and microplates and when performing the test. Protect the Substrate Solution "B" (TMB) from strong light and avoid vibration of the bench surface where the test is undertaken.
6. Upon receipt, store the kit at 2...8°C into a temperature controlled refrigerator or cold room.
7. Do not interchange components between different lots of the kits. It is recommended that components between two kits of the same lot should not be interchanged.
8. Check that the reagents are clear and do not contain visible heavy particles or aggregates. If not, advise the laboratory supervisor to initiate the necessary procedures for kit replacement.
9. Avoid cross-contamination between serum/plasma samples by using disposable tips and changing them after each sample.
10. Avoid cross-contamination between kit reagents by using disposable tips and changing them between the use of each one.
11. Do not use the kit after the expiration date stated on the external container and internal (vials) labels.
12. Treat all specimens as potentially infective. All human serum specimens should be handled at Biosafety Level 2, as recommended by the Center for Disease Control, Atlanta, U.S. in compliance with what reported in the Institutes of Health's publication: "Biosafety in Microbiological and Biomedical Laboratories", ed. 1984. Materials from human origin may have been used in the preparation of the Negative Control of the kit. These materials have been tested with tests kits with accepted

performance and found negative for antibodies to HIV 1/2, HCV, TP and HBsAg. However, there is no analytical method that can assure that infectious agents in the specimens or reagents are completely absent. Therefore, handle reagents and specimens with extreme caution as if capable of transmitting infectious diseases. Bovine derived sera have been used for stabilizing of the positive and negative controls. Bovine serum albumin (BSA) and fetal calf sera (FCS) are derived from animals from BSE/TSE free-geographical areas.

13. The use of disposable plastic-ware is recommended in the preparation of the liquid components or in transferring components into automated workstations, in order to avoid cross contamination.
14. Waste produced during the use of the kit has to be discarded in compliance with national directives and laws concerning laboratory waste of chemical and biological substances. In particular, liquid waste generated from the washing procedure, from residuals of controls and from samples has to be treated as potentially infective material and inactivated before waste. Suggested procedures of inactivation are treatment with a 10% final concentration of household bleach for 16-18 hrs or heat inactivation by autoclave at 121°C for 20 min.
15. Accidental spills from samples and operations have to be adsorbed with paper tissues soaked with household bleach and then with water. Tissues should then be discarded in proper containers designated for laboratory/hospital waste.
16. The Stop Solution contains 0.5M sulphuric acid. Avoid contact with skin and eyes. In the event of contact, rinse immediately with plenty of water.
17. ProClin™ 300 0.1% used as preservative, can cause sensation of the skin. Wipe up spills immediately or wash with water if come into contact with the skin or eyes.
18. Do not smoke, eat, drink or apply cosmetics in areas in which specimens or kit reagents are handled
19. Other waste materials generated from the use of the kit (example: tips used for samples and controls, used microplates) should be handled as potentially infective and disposed according to national directives and laws concerning laboratory wastes.
20. Do not pipette by mouth.

## G. SPECIMEN: PREPARATION AND RECOMMENDATIONS

### 1. Specimen Collection:

No special patient's preparation required. Collect the specimen in accordance with the normal laboratory practice. Either fresh serum or plasma specimens can be used with this assay. Blood collected by venipuncture should be allowed to clot naturally and completely - the serum/plasma must be separated from the clot as early as possible as to avoid haemolysis of the RBC. Care should be taken to ensure that the serum specimens are clear and not contaminated by microorganisms. Any visible particular matters in the specimen should be removed by centrifugation at 3000 RPM (round per minutes) for 20 minutes at room temperature or by filtration.

2. Plasma specimens collected into EDTA, sodium

citrate or heparin may be tested, but **highly lipaemic, icteric, or hemolytic specimens should not be used** as they can give false results in the assay. **Do not heat inactivate specimens.** This can cause deterioration of the target analyte. Samples with visible microbial contamination should never be used.

3. EIAgen HBsAg Kit is intended ONLY for testing of individual serum or plasma samples. Do not use the assay for testing of cadaver samples, saliva, urine or other body fluids, or pooled (mixed) blood.
4. **Transportation and Storage:** Store specimens at 2...8°C. Specimens not required for assaying within 7 days should be stored frozen (-20°C or lower). Multiple freeze-thaw cycles should be avoided. For shipment, samples should be packaged and labeled in accordance with the existing local and international regulations for transportation of clinical samples and ethological agents.

## H. PREPARATION OF COMPONENTS AND WARNINGS

The components of the kit will remain stable through the expiration date indicated on the label and package when stored between 2...8°C, do not freeze. To assure maximum performance of EIAgen HBsAg Kit, during storage, protect the reagents from contamination with microorganism or chemicals.

### 1. Microplates:

Allow the microplate to reach room temperature (18...30°C), about 1 hr, before opening the pouch. Unused strips have to be placed inside the plastic sealable bag, with the desiccant supplied and stored at 2...8°C. After first opening, remaining strips are stable one month at 2...8°C.

### 2. Negative Control:

Ready to use. Mix well on vortex before use.

### 3. Positive Control:

Ready to use. Mix well on vortex before use. Handle this component as potentially infectious.

### 4. Conjugate:

Ready to use. Mix well on vortex before use.

Be careful not to contaminate the liquid with oxidizing chemicals, air-driven dust or microbes. If this component has to be transferred use only plastic, possibly sterile disposable containers.

### 5. Assay Diluent:

Ready to use. Mix well on vortex before use.

### 6. Substrate Solution "A" and "B":

Ready to use. Mix well on vortex before use.

Avoid contamination of the liquid with oxidizing chemicals, air-driven dust or microbes. Do not expose to strong illumination, oxidizing agents and metallic surfaces. If this component has to be transferred use only plastic, and if possible, sterile disposable container.

### 7. Stop Solution:

Ready to use. Mix well on vortex before use.

#### 8. Wash Buffer Concentrate 20x (vial of 30 mL):

The whole content of the 20x concentrated solution has to be diluted with distilled/deionized water up to 600 mL (up to 1000 mL for the vial of 50 mL), the volume is reported on the label, and mixed gently end-over-end before use. As some salt crystals may be present into the vial, take care to dissolve all the content when preparing the solution. During preparation avoid foaming as the presence of bubbles could impact on the efficiency of the washing cycles.

**Note: Once diluted, the wash solution is stable for 1 week at room temperature and two weeks at 2...8°C.**

### I. INSTRUMENTS AND TOOLS USED IN COMBINATION WITH THE KIT

1. Micropipettes have to be calibrated to deliver the correct volume required by the assay and must be submitted to regular decontamination (household alcohol, 10% solution of bleach, hospital grade disinfectants) of those parts that could accidentally come in contact with the sample or the components of the kit. They should also be regularly maintained in order to show a precision of 1% and a trueness of  $\pm 2\%$ . Decontamination of spills or residues of kit components should also be carried out regularly.
2. The ELISA incubator has to be set at  $+37^{\circ}\text{C}$  (tolerance of  $\pm 0.5^{\circ}\text{C}$ ) and regularly checked to ensure the correct temperature is maintained. Both dry incubators and water baths are suitable for the incubations, provided that the instrument is validated for the incubation of ELISA tests.

#### 3. INSTRUCTIONS FOR WASHING

The ELISA washer is extremely important to the overall performances of the assay:

- A good washing procedure is essential in order to obtain correct and precise analytical data.
- It is therefore, recommended to use a good quality ELISA microplate washer, maintained at the best level of washing performances. In general, no less than **5** automatic washing cycles of 350-400  $\mu\text{L}$ /well are sufficient to avoid false positive reactions and high background.
- To avoid cross-contaminations of the plate with specimen or HRP-conjugate, after incubation, do not discard the content of the wells but allow the plate washer to aspirate it automatically.
- Assure that the microplate washer liquid dispensing channels are not blocked or contaminated and sufficient volume of Wash buffer is dispensed each time into the wells.
- In case of manual washing, we suggest to carry out **5 washing cycles**, dispensing 350-400  $\mu\text{L}$ /well and aspirating the liquid for **5** times. If poor results (high background) are observed, increase the washing cycles or soaking time per well.
- In any case, the liquid aspirated out the strips should be treated with a sodium hypochlorite solution at a final concentration of 2.5% for 24 hours, before they are wasted in an appropriate way.
- The concentrated Wash buffer should be diluted **1:20** before use. If less than a whole plate is used, prepare the proportional volume of solution.

4. The ELISA microplate reader has to be equipped with a reading filter of 450nm and ideally with a second filter (630nm) for blanking purposes. Its standard performances should be (a) bandwidth  $\leq 10$  nm; (b) absorbance range from 0 to  $\geq 2.0$ ; (c) linearity to  $\geq 2.0$ ; repeatability  $\geq 1\%$ . Blanking is carried out on the well identified in section "Assay Procedure". The optical system of the reader has to be calibrated regularly to ensure that the correct optical density is measured. It should be regularly maintained according to the manufacturer's instructions.
5. When using an ELISA automated work station, all critical steps (dispensation, incubation, washing, reading, shaking, data handling) have to be carefully set, calibrated, controlled and regularly serviced in order to match the values reported in section O "Internal Quality Control". The assay protocol has to be installed in the operating system of the unit and validated as for the washer and the reader. In addition, the liquid handling part of the station (dispensation and washing) has to be validated and correctly set. Particular attention must be paid to avoid carry over by the needles used for dispensing samples and for washing. This must be studied and controlled to minimize the possibility of contamination of adjacent wells due to strongly reactive samples, leading to false positive results. The use of ELISA automated work stations is recommended for blood screening and when the number of samples to be tested exceed 20-30 units per run. If using fully automated equipment, during incubation, do not cover the plates with the plate cover. The tapping out of the remainders inside the plate after washing, can also be omitted.

### L. PRE ASSAY CONTROLS AND OPERATIONS

1. Check the expiration date of the kit printed on the external label of the kit box. Do not use if expired.
2. Check that the liquid components are not contaminated by naked-eye visible particles or aggregates. Check that the Substrate Solutions "A" and "B" are colorless or pale blue by aspirating a small volume of it with a sterile transparent plastic pipette. Check that no breakage occurred in transportation and no spillage of liquid is present inside the box. Check that the aluminum pouch, containing the microplate, is not punctured or damaged.
3. Dilute all the content of the 20x concentrated Wash Buffer as described above.
4. Allow all the other components to reach room temperature ( $18...30^{\circ}\text{C}$ ), about 1 hr and then mix as described.
5. Set the ELISA incubator at  $+37^{\circ}\text{C}$  and prepare the ELISA washer by priming with the diluted wash solution, according to the manufacturers instructions. Set the right number of washing cycles as found in section I.3.
6. Check that the ELISA reader has been turned on at least 20 minutes before reading.
7. If using an automated workstation, turn it on, check settings and be sure to use the right assay protocol.
8. Check that the micropipettes are set to the required volume.

9. Check that all the other equipment is available and ready to use.
10. In case of problems, do not proceed further with the test and advise the supervisor.

### M. ASSAY PROCEDURE

The assay has to be carried out according to what reported below, taking care to maintain the same incubation time for all the samples in testing.

#### MANUAL ASSAY:

1. **Preparation:** Mark three wells as Negative control (e.g. B1, C1, D1), two wells as Positive control (e.g. E1, F1) and one Blank (e.g. A1, neither samples nor Conjugate should be added into the Blank well). If the results will be determined by using dual wavelength plate reader, the requirement for use of Blank well could be omitted. Use only number of strips required for the test.
2. **Adding Diluent:** Add 20µL of Assay Diluent into each well except the Blank.
3. **Adding Sample:** Add 100µL of Positive control, Negative control, and Specimen into their respective wells except the Blank. **Note: Use a separate disposal pipette tip for each specimen, Negative Control, Positive Control to avoid cross-contamination. Mix by tapping the plate gently.**
4. **Incubating:** Cover the plate with the plate cover and incubate for **60 minutes at 37°C**.
5. **Adding Conjugate:** At the end of the incubation, remove and discard the plate cover. Add 50µL Conjugate into each well except the Blank, and mix by tapping the plate gently.
6. **Incubating:** Cover the plate with the plate cover and incubate for **30 minutes at 37°C**.
7. **Washing:** At the end of the incubation, remove and discard the plate cover. Wash each well **5 times** with diluted Washing buffer. Each time allow the microwells to soak for **30-60 seconds**. After the final washing cycle, turn down the plate onto blotting paper or clean towel and tap it to remove any remainders (see section I.3).
8. **Coloring:** Add 50µL of Substrate Solution "A" and 50µL Substrate Solution "B" into each well including the Blank. Incubate the plate at **37°C for 30 minutes avoiding light**. The enzymatic reaction between the Substrate Solutions and the Conjugate produces blue color in Positive control and HBsAg positive sample wells.
9. **Stopping Reaction:** Using a multichannel pipette or manually, add 50µL Stop solution into each well and mix gently. Intensive yellow color develops in Positive control and HBsAg positive sample wells.
10. **Measuring the Absorbance:** Calibrate the plate reader with the Blank well and read the absorbance at **450nm**. If a dual filter instrument is used, set the reference wavelength at **630nm**. Calculate the Cut-off value and evaluate the results. (**Note:** read the absorbance within **10 minutes** after stopping the reaction).

#### Important notes:

1. If the second filter is not available, ensure that no finger prints are present on the bottom of the

microwell before reading at 450nm. Finger prints could generate false positive results on reading.

2. Reading has should ideally be performed immediately after the addition of the Stop Solution but definitely no longer than 10 minutes afterwards. Some self oxidation of the substrate can occur leading to a higher background.

### N. ASSAY SCHEME

Steps	Operations
Assay Diluent	20 µL
Controls and Samples	100 µL
<b>1<sup>st</sup> incubation</b>	<b>60 min</b>
Temperature	+37°C
Conjugate	50 µL
<b>2<sup>nd</sup> incubation</b>	<b>30 min</b>
Temperature	+37°C
Wash Step	5 cycles (see section I.3)
Substrate "A"	50 µL
Substrate "B"	50 µL
<b>3<sup>rd</sup> incubation</b>	<b>30 min (avoiding light)</b>
Temperature	+37°C
Stop Solution	50 µL
Reading OD	450/630nm

An example of dispensation scheme is reported below (valid for both incubation time procedures):

		Microplate											
		1	2	3	4	5	6	7	8	9	10	11	12
A	BLK	S3											
B	NC	S4											
C	NC	S5											
D	NC	S6											
E	PC	S7											
F	PC	S8											
G	S1	S9											
H	S2	S10											

Legenda: BLK = Blank NC = Negative Control  
PC = Positive Control S = Sample

### O. INTERNAL QUALITY CONTROL

The test results are valid if the Quality Control criteria are fulfilled. It is recommended that each laboratory must establish appropriate quality control system with quality control material similar to or identical with the patient sample being analyzed.

- The A value of the Blank well, which contains only Chromogen and Stop Solution, is < 0.080 at 450 nm.
- The A values of the Positive Control must be ≥ 0.800 at 450/630nm or at 450nm after blanking.
- The A values of the Negative Control must be < 0.100 at 450/630nm or at 450nm after blanking.

If one of the Negative control A values does not meet the Quality Control criteria, it should be discarded and the mean value calculated again using the remaining two values. If more than one Negative control A values do not meet the Quality Control Range specifications, the test is invalid and must be repeated.

If the results of the test match the requirements stated above, proceed to the next section. If they do not, do not proceed any further and perform the following checks:

Problems	Check
<b>Blank well</b> ≥ 0.080 OD at 450nm	1. that the Substrate Solutions have not become contaminated during the assay
<b>Negative Control (NC)</b> ≥ 0.100 OD at 450/630nm or at 450nm after blanking	1. that the washing procedure and the washer settings are as validated in the pre qualification study; 2. that the proper washing solution has been used and the washer has been primed with it before use; 3. that no mistake has been done in the assay procedure (dispensation of positive control instead of the negative one); 4. that no contamination of the negative control or of the wells where the control was dispensed has occurred due to spills of positive samples or of the conjugate; 5. that micropipettes have not become contaminated with positive samples or with the conjugate 6. that the washer needles are not blocked or partially obstructed.
<b>Positive Control</b> < 0.800 OD at 450/630nm or at 450nm after blanking	1. that the procedure has been correctly performed; 2. that no mistake has occurred during the distribution of the control (dispensation of negative control instead of positive control); 3. that the washing procedure and the washer settings are as validated in the pre qualification study; 4. that no external contamination of the positive control has occurred.

If any of the above problems have occurred, report the problem to the supervisor for further actions.

### P. RESULTS

Each microplate should be considered separately when calculating and interpreting the results of the assay, regardless of the number of plates concurrently processed. The results are calculated by relating each specimen absorbance (A) value to the Cut-off value (C.O.) of the plate. If the Cut-off reading is based on single filter plate reader, the results should be calculated by subtracting the Blank well A value from the print report values of specimens and controls. In case the reading is based on dual filter plate reader, do not subtract the Blank well A value from the print report values of specimens and controls.

$$\text{Cut-Off (C.O.)} = \text{NC mean} + 0.06$$

The value found for the test is used for the interpretation of results as described in the next paragraph.

An example of calculation is reported below:

### Example:

#### 1. Quality Control

Blank well A value: A1= 0.025 at 450nm (Note: blanking is required only when reading with single filter at 450nm)

**Well No.:** B1 C1 D1  
 Negative control A values after blanking: 0.020 0.012 0.016

**Well No.:** E1 F1  
 Positive control A values after blanking: 2.421 2.369  
 All control values are within the stated quality control range

$$2. \text{ Calculation of Nc: } = \frac{(0.020+0.012+0.016)}{3} = 0.016$$

$$3. \text{ Calculation of the Cut-off: (C.O.)} = 0.016 + 0.06 = 0.076$$

### Q. INTERPRETATION OF RESULTS

#### Negative Results (A / C.O. < 0.9):

Specimens giving absorbance less than the Cut-off value are negative for this assay, which indicates that no hepatitis B virus surface antigen has been detected with EIAGEN HBsAg Kit, therefore the patient is probably not infected with HBV and the blood unit do not contain hepatitis B virus surface antigen and could be transfused in case that other infectious diseases markers are also absent.

#### Positive Results (A / C.O. > 1.1):


Specimens giving an absorbance equal to or greater than the Cut-off value are considered initially reactive, which indicates that hepatitis B virus surface antigen has probably been detected using EIAGEN HBsAg Kit. All initially reactive specimens should be retested in duplicates using EIAGEN HBsAg Kit before the final assay results interpretation. Repeatedly reactive specimens can be considered positive for hepatitis B virus surface antigen with EIAGEN HBsAg Kit.

#### Borderline (A / C.O. = 0.9-1.1):

Specimens with absorbance to Cut-off ratio between 0.9 and 1.1 are considered borderline and retesting of these specimens in duplicates is required to confirm the initial results.

**Follow-up, confirmation and supplementary testing of any positive specimen with other analytical system (e.g. PCR) is required. Clinical diagnosis should not be established based on a single test result. It should integrate clinical and other laboratory data and findings.**

#### INITIAL RESULTS INTERPRETATION AND FOLLOW-UP ALL INITIALLY REACTIVE OR BORDERLINE SAMPLES

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If, after retesting of the initially reactive samples, both wells are negative results (A/C.O.<0.9), these samples should be considered as non-repeatable positive (or false positive) and recorded as negative. As with many very sensitive ELISA assays, false positive results can occur due to the several reasons, most of which are connected with, but not limited to, inadequate washing step. For more information regarding "Troubleshooting Guide" (section S).

If after retesting in duplicates, one or both wells are positive results, the final result from this ELISA test should be recorded as repeatedly reactive. Repeatedly reactive specimens could be considered positive for hepatitis B virus surface antigen and therefore the patient is probably infected with HBV and the blood unit must be discarded.

After retesting in duplicates, samples with values close to the Cut-off value should be interpreted with caution and considered as "borderline" zone sample, or uninterpretable for the time of testing.

**Important notes:**

1. Interpretation of results should be done under the supervision of the laboratory supervisor to reduce the risk of judgment errors and misinterpretations.
2. When test results are transmitted from the laboratory to another department, attention must be paid to avoid erroneous data transfer.
3. Diagnosis of viral hepatitis infection has to be taken and released to the patient by a suitably qualified medical doctor.

**R. PERFORMANCES**

Evaluation of Performances has been conducted in accordance to what reported in the Common Technical Specifications or CTS (art. 5, Chapter 3 of IVD Directive 98/79/EC).

Evaluation studies carried out in Paul-Ehrlich-Institut (PEI), German Red Cross Institute Baden-Württemberg – Hessen, and three blood banks, demonstrated the following performance characteristics of EIAGEN HBsAg Kit.

**1. SPECIFICITY**

When evaluated on European blood donors (n=5038), the overall diagnostic specificity of the kit was 99.78%.

During multi-center evaluation (Site A, B and C), EIAGEN HBsAg Kit demonstrated specificity of 99.92%.

Laboratory	Number	EIAGEN HBsAg Kit		
		-	+	Specificity
"A" blood bank	1958	1955	3	99.85%
"B" blood bank	2518	2516	2	99.92%
"C" blood bank	6344	6340	4	99.94%
Total	10820	10811	9	99.92%

**2. SENSITIVITY**

EIAGEN HBsAg Kit was evaluated for sensitivity on 22 HBV commercial available HBV seroconversion panels, and on total 403 HBsAg positive including 146 HBsAg HBV genotyped and HBsAg subtyped plasma samples available at the Paul-Ehrlich-Institut. With respect to seroconversion sensitivity, the results for EIAGEN HBsAg Kit on the 22 HBV seroconversion panels showed a sensitivity level at least equivalent with the range of current CE marked HBsAg screening assays for which PEI holds data. 10 additional seroconversion panels were tested in-house. The seroconversion sensitivity was comparable to other CE-marked HBsAg screening test. With respect to diagnostic sensitivity EIAGEN HBsAg Kit detected all positive samples as positive, including the HBV genotypes A-F or HBsAg subtypes examined.

In conclusion, the overall score of EIAGEN HBsAg Kit for the seroconversion sensitivity was comparable with other CE marked HBsAg test kits for which PEI holds data and all 403 HBsAg positive samples were reactive giving an overall sensitivity of 100%.

**3. ANALYTICAL SENSITIVITY**

0.1 IU/mL (NIBSC 00/588)

**4. ANALYTICAL SPECIFICITY**

No interference was observed with samples from patients with high-level of rheumatoid factor, and pregnant woman. Same day and frozen specimens have been tested to check for interferences due to collection and storage. Total of 100 samples reactive for anti-HBc, anti-HCV and anti-HIV-1 were screened for HBsAg with EIAGEN HBsAg Kit. 98 out of 100 samples were negative for HBsAg. 200 blood samples from patients were also tested with EIAGEN HBsAg Kit. 191 out of 200 samples had negative screening results for HBsAg. 8 out of 9 samples with initial reactive screening results had repeat reactive test results with EIAGEN HBsAg Kit but hepatitis B virus was not confirmed in all cases.

**5. DETECTION OF MUTATIONS**

Panel of 108 samples sequenced by PCR were tested to demonstrate the performance of EIAGEN HBsAg Kit in detection of HBsAg mutations. The results are given in the table below.

Background	Number	EIAGEN HBsAg Kit
adr (+)	wild type	35
	4 mutations	5
adw (+)	wild type	37
	16 mutations	25
ayw (+)	wild type	2
	2 mutations	2
ayr (+)	2 mutations	2
Total	108	101

## 6. REPRODUCIBILITY

Reproducibility		Within Run		Between Run	
Specimen Type	n	Mean OD	CV%	Mean OD	CV%
0.1 IU/mL HBsAg	10	0.155	10.6%	0.150	11.0%
Weak positive	10	0.457	9.0%	0.432	9.5%
Moderate positive	10	1.572	7.0%	1.437	7.5%
Strong positive	10	2.327	4.2%	2.302	4.4%
Positive control	10	2.322	4.1%	2.315	4.2%

## S. SUGGESTIONS FOR TROUBLESHOOTING

Adherence to assay procedure and specifications, as well as a correct use of reagents and proper pipetting, may help to avoid the following kinds of errors:

ERROR	POSSIBLE CAUSES / SUGGESTIONS
OD very different ( $\pm 50\%$ ) from OD reported on QC	<ul style="list-style-type: none"> <li>- incorrect dispensing volume of reagents (suggestion: check the correspondence between the volume dispensed by the pipette and the one required by the assay; re-calibrate again pipettes)</li> <li>-incorrect temperature or incorrect incubation time (suggestion: more care in the incubator maintenance; note down the beginning of the incubation)</li> <li>-error in washing or in photometer reading (suggestion: check operating or settings of respective instruments)</li> <li>-contamination of Substrate Solutions or Conjugate (suggestion: use only disposable and clean plastic containers)</li> </ul>
Low reproducible results	<ul style="list-style-type: none"> <li>-not constant dispensing volume of samples or reagents (suggestion: check the pipettes precision and the correspondence between the volume dispensed by the pipette and the one required by the assay; re-calibrate again pipettes)</li> <li>-error in washing or in reading (suggestion: check operating or settings of respective instruments)</li> <li>-contamination of Substrate Solutions (suggestion: use only disposable and clean plastic containers)</li> <li>-pollution or degradation of reagents (suggestion: use appropriate tips, disposable and clean plastic containers for reagents and high quality distilled or equivalent water)</li> </ul>
no colorimetric reaction after addition of substrate solutions	<ul style="list-style-type: none"> <li>-some reagent not pipetted</li> <li>- strong contamination of Conjugate or Substrate Solutions</li> <li>-errors in performing the assay procedure (e.g. accidental pipetting of reagents in a wrong sequence or from the wrong vial, etc.)</li> </ul>
too low reaction (too low ODs)	<ul style="list-style-type: none"> <li>-incubation time too short, incubation temperature too low</li> <li>-incorrect conjugate dilution</li> </ul>
too high reaction (too high ODs)	<ul style="list-style-type: none"> <li>-incorrect conjugate dilution</li> <li>-incubation time too long, incubation temperature too high</li> <li>-water quality for wash buffer insufficient (low grade of deionization)</li> <li>-insufficient washing (conjugates not properly removed)</li> </ul>
unexplainable outliers	<ul style="list-style-type: none"> <li>-contamination of pipettes, tips or containers</li> <li>-inconstant and insufficient washing (conjugates not properly removed)</li> </ul>

too high within-run CV%	<ul style="list-style-type: none"> <li>-reagents and/or strips not pre-warmed to Room Temperature prior to use</li> <li>- plate washer is not washing correctly (suggestion: clean washer head)</li> </ul>
too high between-run CV%	<ul style="list-style-type: none"> <li>-incubation conditions not constant (time, temperature)</li> <li>-controls and samples not dispensed at the same time (with the same intervals) (check pipetting order)</li> <li>-person-related variation</li> </ul>

## T. AUTOMATION

The procedures identified in this Instruction for Use are for manual testing only. When using automated instruments, follow the procedures that are contained in the operator's manual provided by the device manufacturer. Laboratories must follow their approved validation procedures to demonstrate compatibility of this product on automated systems.

## U. LIMITATIONS

1. Positive results must be confirmed with another available method and interpreted in conjunction with the patient clinical information.
2. Antigens may be undetectable during the early stage of the disease. Therefore, negative results obtained with EIAgen HBsAg Kit are only indication that the sample does not contain detectable level of hepatitis B virus surface antigen and any negative result should not be considered as conclusive evidence that the individual is not infected with HBV or the blood unit is not infected with HBV.
3. If, after retesting of the initially reactive samples, the assay results are negative, these samples should be considered as non-repeatable (false positive) and interpreted as negative. As with many very sensitive ELISA assays, false positive results can occur due to the several reasons, most of which are related but not limited to inadequate washing step. For more information please refer to "Troubleshooting Guide", or contact Adaltis technical support for further assistance.
4. The most common assay mistakes are: using kits beyond the expiry date, bad washing procedures, contaminated reagents, incorrect assay procedure steps, insufficient aspiration during washing, failure to add specimens or reagents, improper operation with the laboratory equipment, timing errors, the use of highly hemolyzed specimens or specimens containing fibrin, incompletely clotted serum specimens.
5. The prevalence of the marker will affect the assay's predictive values.
6. This assay cannot be utilized to test pooled (mixed) plasma. EIAgen HBsAg Kit has been evaluated only with individual serum or plasma specimens.
7. EIAgen HBsAg Kit is a qualitative assay and the results cannot be used to measure antigen concentration.

## 8. INDICATIONS OF INSTABILITY DETERIORATION OF THE REAGENT:

Values of the Positive or Negative controls, which are out of the indicated quality control range, are indicators of possible deterioration of the reagents and/or operator or equipment errors. In such case, the results should be considered as invalid and the


samples must be retested. In case of constant erroneous results and proven deterioration or instability of the reagents, immediately substitute the reagents with new one or contact Adaltis technical support for further assistance.

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